

2. Performance characteristics of tests used in the initial evaluation of patients at risk of renal disease

Date written: February 2003

Date submitted: June 2004

GUIDELINES

No recommendations possible

BACKGROUND

End-stage renal disease (ESRD) develops in about 1500 people each year in Australia, of whom about 1000 are over 50 years of age (a risk of about 200 per million per year). Diabetes, hypertension, impaired renal function and proteinuria are likely independent risk factors for the development of ESRD. About 5% of the general population have proteinuria and these individuals are approximately 15 times more likely to develop ESRD than those without proteinuria. Because of the absence of well done cohort studies, estimating an individual's risk of developing ESRD using values for proteinuria, blood pressure, renal function and diabetes is not possible. The aim of this guideline is to review the different tests available for assessing proteinuria and albuminuria.

SEARCH STRATEGY

Databases searched: MeSH terms (when available) and text words for each intervention were combined with text words (no MeSH terms exist) for the comparison. These were further combined with MeSH terms and text words for the outcomes. The search was carried out in Medline (1966 – January Week 1 2003).
Date of search: 21 January 2003.

WHAT IS THE EVIDENCE?

There are no randomised controlled trials on this topic.

SUMMARY OF THE EVIDENCE

Not possible.

SUGGESTIONS FOR CLINICAL CARE

(Suggestions are based on level III and IV sources)

- **Dipstick protein has been identified to have poor sensitivity and specificity. Dipstick albumin is more reliable but operator (and technique) dependant.**

- For initial testing for albuminuria, first morning urine albumin sample is preferred, however, a random urine albumin creatinine ratio (ACR) is acceptable. If abnormal, one or more timed urines are required to confirm microalbuminuria.

- Protein creatinine ratio (PCR) is an accurate test for diagnosis of significant proteinuria (>1 g/day) but timed urines are required to establish a baseline.

Dipstick tests

Dipstick testing for proteinuria

The performance characteristics of dipstick testing have been compared to the reference standard of 24 h urine protein quantification. Six relevant studies containing 1558 patients were identified and assessed by meta-analysis by Barratt *et al.* (see Table 1).¹ The studies were of variable quality and carried out in high prevalence populations (mostly renal, rheumatology and obstetric outpatient settings) but use similar tests (mostly Multistix) and similar reference and test thresholds (300 mg and 1+, respectively). Using this test to detect proteinuria, there was a trade off between sensitivity and specificity:

Sensitivity: 80% Specificity: 85%

Sensitivity: 90% Specificity: 67%

Sensitivity: 60% Specificity: 95%.

The populations used in these studies had high prevalences of proteinuria, mostly over 30% (although the largest had only 7%). It is unlikely that dipstick testing will be employed to screen populations with a prevalence of >30% beyond the first round of screening. Hence, using dipsticks as a screening tool for proteinuria has marked limitations. These include:

- poor specificity and sensitivity yielding high false positive and false negative rates
- the generally poor quality of studies in this field and absence of studies in suitable at-risk populations, and
- the likely incidence of proteinuria in at-risk populations of <10%.

Dipstick testing for microalbuminuria

Studies found using the search criteria that related to spot urine testing for microalbuminuria and used 24 h urine collection for albuminuria as gold standard were included. Abstracts of all relevant studies were reviewed and seven articles were identified, with two of these being excluded (one due to poor quality and one due to failure to use appropriate reference standard). Of the five studies reviewed (see Table 2), the overall standard was poor, with the exception of Gilbert *et al.*² All studies tested the Micral-Test (Boehringer Mannheim, GmbH Mannheim, Germany) test strips and only one used non-diabetic subjects (hypertensive patients).³

Other studies suggest that the Micral-Test is operator dependent with superior results seen in some professional groups and after training in the test technique. The test strips must be immersed in the urine for 5 s and read 5 min later. Other studies have suggested that test accuracy is very dependent on rigorous adherence to the timing aspects of the test.

It appears that the Micral-Test has better sensitivity and specificity for the detection of microalbuminuria than the dipstick for proteinuria. It has also been studied in appropriate populations, although the quality of these studies still leaves a lot to be desired. The prevalence in the study populations is very high and would not be the case beyond the first round of screening, making the sensitivity and specificity more important.

General comments about dipstick testing

The studies all suggest that using higher cut-offs than 1+ (for proteinuria) and 20 mg/L (for microalbuminuria) would improve the sensitivity and specificity, at the risk of missing those with very low grade proteinuria or albuminuria. The outcome studies all suggest that those patients with higher levels of proteinuria are at higher risk of renal deterioration and are more likely to benefit from treatment. Therefore, using higher cut-offs for detection may make the urine tests more viable, allowing clinicians to take advantage of the immediacy, ease and cost savings of dipstick testing with a low risk of not treating those patients most likely to benefit from the interventions planned.

Random urine testing for protein and albumin

The performance characteristics of random urine testing have been compared to the reference standard of 24 h urine collections. The quality of the evidence is marred by methodological and analytical flaws. Assessment of adequacy of 24 h collections via determination of the creatinine excretion is variable. In studies where the adequacy of collection is assessed, up 15%–23% of samples were excluded from analysis.^{4,5} The patient populations

Table 1 Studies that compared dipstick to 24 h urine protein¹

Author, Year	Setting	Reference threshold	Test type	Test threshold	N	Sensitivity (%)	Specificity (%)	Prevalence (%)	PPV	NPV
Ralston, 1988	Rheumatology outpatient department	300 mg	Multistix	1+	102	100	36	42	47	100
Harrison, 1989	Renal admissions	300 mg	N-Multistix-SG	1+	100	82	77	49	77	82
Allen, 1991	Sequential inpatients & outpatients	250 mg	Multistix-10SG	1+	136	47	97	46	93	68
Meyer, 1994	Hypertensive pregnancies	300 mg	Multistix	1+	300	67	74	81	92	34
Brown, 1995	Hypertensive pregnant inpatients	300 mg	Multistix-10SG	1+	230	86	39	30	38	87
Highy, 1995	Pregnant women	300 mg	Micro-bumintest	1+	690	69	86	7	72	98

Table 2 Studies that compared dipstick to 24 h urine albumin

Author, Year	Reference threshold	Test type	Test threshold	No.	Sensitivity (%)	Specificity (%)	Prevalence (%)	PPV	NPV
Gilbert, ² 1992	20 mg/L	Micral	20 mg/L	298	92	92	33	85	96
Gerber, ³ 1998	20 mg/L	Micral	20 mg/L	171	81	89	22	68	94
Spooren, ²¹ 1992	20 µg/min	Micral	20 mg/L	132	82?	86?	?		
Marshall, ²² 1992	20 mg/L	Micral	20 mg/L	112	100	91	30	83	100
Mogensen, ²³ 1997	20 mg/L	Micral II	20 mg/L	2228	97	71	52	78	96

are mostly taken from specialist hospital clinics and involve both inpatients and outpatients, and may include multiple samples from the same patient.

Random urine testing for protein

Urine protein concentration: protein creatinine ratio

Good correlation between PCR and 24 h urine total protein (UTP) has been established in a wide range of patient groups (e.g. normal, known renal disease, diabetes). A meta-analysis of 7 studies (marked*) containing 873 patients found the sensitivity to be 95% (92%–97%) and specificity 91% (87%–94%)¹ (see Table 3).

Protein-osmolality ratio has also been used to correct for urine concentration with similar degrees of sensitivity (91%–96%) and specificity (93%–98%).⁶

Random urine testing for albumin

Review of the literature is made difficult by the various definitions of microalbuminuria, wide variations in test thresholds, and lack of standard comparator. Some studies compared the samples from the same collection, others from samples taken weeks apart. Some studies used overnight urine albumin excretion (which has also been shown to predict development of proteinuria in diabetics). There is some difficulty in comparing overnight samples to random daytime samples as the relationship between daytime and overnight albumin excretion is variable. Overnight albumin excretion is 25% lower than the 24 h albumin excretion rate (AER).⁷ The intra-subject day-to-day variability in albumin excretion in a daytime sample is also greater than in an overnight sample (presumably due to the effect of variable exercise). Therefore, studies using 24 h urine collections are the gold standard but other studies using timed overnight urines are included (see Table 4).

Urine albumin concentration (UAC)

The best correlation with AER is found with a first morning sample with a sensitivity from 86%–97% and specificity from 74%–97%. Random urine albumin concentration correlates poorly with AER in a timed specimen and lacks sensitivity (56%–96%) and specificity (81%–80%).

Albumin: creatinine ratio (ACR)

Correction for urine flow rate increases the sensitivity and specificity with best results seen with early morning samples. However, almost all of the studies use the same sample for assessment of ACR and timed urine and this is a serious methodological flaw.

A systematic review of the literature in urine albumin testing in diabetics has been performed by Scheid *et al.*⁸ Ten studies comparing ACR to timed albumin collection and containing 1691 patients were reviewed. Due to the

heterogeneity of studies, the trials were not pooled and data is presented as the range (see below).

	No.	Sensitivity	Specificity
First morning ACR	355	70–97	91–94
Morning ACR	240	71–100	95–98
Random ACR	1096	56–97	81–92

Limitations of protein creatinine ratio/albumin creatinine ratio

The use of urine creatinine excretion to correct for variation in urine concentration needs to be applied with an understanding of its limitations. Gender, age and renal function will all affect the creatinine excretion. Unless this is considered, there is a tendency to underestimate proteinuria in men and overestimate its prevalence in the elderly of both genders due to differences in the amount of creatinine excreted with different muscle mass.^{9,10} This difference has significance for screening populations. The NHANES III study found a similar prevalence of microalbuminuria (MA), when defined as UAC > 30 mg/L in men and women (10.9% and 11.8%, respectively) but when applying the same reference range in both genders, found that MA (ACR > 30 mg/g) was more common in women (9.7%) compared with men (6.1%).^{11,12} Gender-specific ranges have been studied in 218 healthy controls and values corresponding to 95% of the respective gender distributions derived. The result is lower in men – 17 mg/mg (1.9 mg/mmol) compared with 25 mg/mg (2.9 mg/mmol) in women. When these different ranges were applied to the NHANES III cohort, the prevalence of MA was the same in the genders.¹³

The excretion of protein and albumin is not constant over 24 h – it is impacted by posture, activity, protein intake and haemodynamic factors. There is evidence that urine protein excretion has a circadian pattern, being highest in the afternoon. There is little corresponding variation in creatinine excretion. Therefore the timing of the sample will markedly influence the predicted proteinuria (mean variation 42%). The best estimates are probably on morning samples (not first void).^{14–16} The variability is less prominent in patients with massive proteinuria. The variability of albumin excretion is similar in magnitude.¹⁷ The average intraindividual daily UAE is 40% and is most marked at low levels of albumin excretion (<30 mg/day).

There is day-to-day variability of excretion of protein and albumin. Koopman *et al.* repeated 24 h urine protein measurements on three consecutive days and found that the variability, expressed as a coefficient of variation (CV), was 15% and not more reliable when assessed at any particular time point.¹⁶ Cundy *et al.* measured 24 h UAE and ACR in 30 diabetics and repeated the test within a month – the CV for 24 h UAE was 98% and the ACR CV was 90% – with poorer agreement at low levels

Table 3 Characteristics of protein creatinine ratio compared to 24 h urine protein

Author	Year	Characteristics	No.	Reference threshold	Test threshold	Sensitivity (%)	Specificity (%)	Correlation
Ginsberg ^{14*}	1983	Renal outpatient department	46	1 g/1.73 m ²	1	93	94	0.97
Parag ^{24*}	1986	Renal outpatient department	34	3000 mg	3000	91	82	0.90
Ralston ^{25*}	1988	Rheumatology inpatient + outpatient department	102	300 mg	0.04	97	97	0.92
Wilson ^{6*}	1993	Normal	131	156 mg men, 90 mg women	0.05	96	90	0.91
Mitchell ^{15*}	1993	Renal outpatient department	52	1000 mg	1000	100	94	0.98
Robert ^{26*}	1997	>65 years, renal outpatient department	71	300 mg	?	93	90	0.64
Ruggenenti ^{27*}	1998	Hypertension/pregnant Trial	98	368 mg		94	92	0.63
Ramos ⁷	1999	Pregnant/hypertensive	47	300 mg	0.5	96	96	0.94
Zelmanovitz ²⁸	1998	Diabetics	167	541 mg	UPC [†]	100	92.9	
					PCR [‡]	100	76.2	

Note: Test thresholds are calculated in different ways in different studies; [†]UPC = urine protein concentration; [‡]PCR = protein creatinine ratio.

Table 4 Characteristics of urinary albumin testing compared to 24 h urine protein

Author, Year	No.	Patient group	Comparator	Reference threshold	Test threshold	Sensitivity (%)	Specificity (%)	Correlation
Random UAC								
Gatling, ²⁹ 1985	159	Diabetic in GP	Overnight	30 µg/min	mg/L 25	56	81	0.45
Watts, ³⁰ 1989	160	Type I diabetics	Overnight	15 µg/min	15	96	80	
				30 µg/min	26	100	85	
Zelmanovitz, ³¹ 1997	95	Diabetic OPD*	24 h	20 µg/min	16.9	100	79	0.91
Ahn, ³² 1999	105	Diabetic OPD	24 h	20 µg/min	31	77	82	0.81
Early morning UAC								
Gatling, ²⁹ 1985	175	Diabetic in GP	Overnight	30 µg/min	mg/L 20	86	97	0.9
Hutchison, ³³ 1988	261	Diabetic hospital OPD	Overnight	30 µg/min	17	96.8	90.9	0.90
Wiegmann, ³⁴ 1990	135	Diabetics	24 h	20 µg/min	30	71	23	
Bakker, ³⁵ 1999	2394	Diabetic	Overnight	20 µg/min	15	88.8	90.2	
Random ACR								
Gatling, ²⁹ 1985	463	Diabetic in GP	Overnight	30 µg/min	mg/mmol 3	80	81	0.43
Nathan, ³⁶ 1987	35	Diabetic OPD	24 h urine	15 µg/min	3.4	94	96	0.82
Zelmanovitz, ³¹ 1997	95	Diabetic OPD	24 h	20 µg/min	1.72	100	74	0.92
Wiegmann, ³⁴ 1990	135	Diabetics	24 h	20 µg/min	3.4	82	81	
Watts, ³⁰ 1989	160	Type I diabetics	Overnight	15 µg/min	2.9	100	80	
				30 µg/min	8	100	96	
Ahn, ³² 1999	105	Diabetic OPD	24 h	20 µg/min	3.7	77	92	
James, ³⁷ 1995	33	>55 years hypertensive OPD non-diabetic	24 h	30 mg/day	3	92	90	
Early morning ACR								
Hutchison, ³³ 1988	261	Diabetic hospital OPD	Overnight	30 µg/min	3	96.8	93.9	0.92
Gatling, ²⁹ 1985	171	Diabetic OPD	Overnight	30 µg/min	3.5	100	95	0.9
Marshall, ³⁸ 1986		Diabetics		30 µg/min	3.5	98	69	
Eshoj, ³⁹ 1987	54	Diabetics	24 h	20 µg/min	3.5	76	97	
Bakker, ³⁵ 1999	2394	Diabetic	Overnight	20 µg/min	2.5	98.5	86.7	

OPD = outpatient department

Drawing conclusions from these studies is difficult given the variability in methodology and thresholds, but it would appear that ACR (preferably a first morning specimen) can provide a useful semi-quantitative measure which is useful in screening for albuminuria.

of albuminuria for both tests.¹⁸ Although the limits of agreement were wide between methods, the correlation was very strong for all comparisons ($R = 0.94\text{--}0.97$). The authors concluded that the variability in albumin excretion was so great that it overcame the poor agreement between tests.

False-negative results for microalbuminuria may occur in the setting of heavy proteinuria. This results from the prozone or Hook effect, which is an inherent problem with antibody-based immunoassays when the antigen (albumin) is in excess.¹⁹ It is standard laboratory practice and recommended by the product inserts for the various assay kits, that urine protein be assessed first, usually by dipstick and samples with heavy proteinuria be diluted to achieve appropriate concentration prior to the assay.

Recent studies have demonstrated that urinary albumin is composed of both immunoreactive and immuno-unreactive fractions.²⁰ When albuminuria is assessed by radio-immuno assay (RIA), as is standard practice, the immuno-unreactive component will be missed. The exact nature of immuno-unreactive albumin is unknown but may involve a conformational change which is not recognized by the antibodies used in the conventional RIA. Total urine albumin can be measured by high performance liquid chromatography (HPLC). No difference is detected in normal controls, however, in diabetics, 91% of patients have significantly more albuminuria when assessed by HPLC. The difference is greatest at low level albuminuria, with a two-fold increase in the number of patients in the microalbuminuria range when assessed by HPLC compared with RIA. The HPLC assay also allows earlier detection of patients who progress from normoalbuminuria to microalbuminuria.²⁰ The assay is not widely available and further studies are required to confirm these findings. In addition, the impact of various interventions such as blockade of the renin-angiotensin system on the immuno-unreactive fraction are unknown.

WHAT DO THE OTHER GUIDELINES SAY?

Kidney Disease Outcomes Quality Initiative:

Under most circumstances, untimed ('spot') urine samples should be used to detect and monitor proteinuria in adults and children.

It is usually not necessary to obtain a timed urine collection (overnight or 24-h) for these evaluations in either children or adults.

First morning specimens are preferred, but random specimens are acceptable if first morning specimens are not available.

In most cases, screening with urine dipsticks is acceptable for detecting proteinuria.

Standard urine dipsticks are acceptable for detecting increased total urine protein.

Albumin-specific dipsticks are acceptable for detecting albuminuria.

Patients with a positive dipstick test (1+ or greater) should undergo confirmation of proteinuria by a quantitative measurement (protein-to-creatinine ratio or albumin-to-creatinine ratio) within three months.

Patients with two or more positive tests temporally spaced by 1–2 weeks should be diagnosed as having persistent proteinuria and undergo further evaluation and management for chronic kidney disease.

Specific guidelines for adults:

When screening adults at increased risk for chronic kidney disease, albumin should be measured in a spot urine sample using either:

Albumin-specific dipstick

Albumin-to-creatinine ratio.

British Diabetes Association:

Everyone with diabetes aged >11 years should be screened annually:

- Screen when free of acute intercurrent illness
- Collect an early morning urine sample
- Use a sensitive assay, specific for albumin
- Measure urine albumin concentration or albumin:creatinine ratio and serum creatinine

Classify as:

Lower risk (absence of microalbuminuria and proteinuria and normal serum creatinine)

Higher risk (microalbuminuria or proteinuria or elevated serum creatinine)

- If lower risk screen annually
- If higher risk, obtain two further urine tests for confirmation, within three months.

Australian Diabetes Society Position Statement (1993):

Testing for microalbuminuria should be performed at least yearly in all diabetic patients.

The finding of a single albumin excretion rate (AER) value in the rate of 20–200 $\mu\text{g}/\text{min}$ is an indication for two additional measurements of AER in the next 6 weeks.

Two of three AER measurements in the range of 20–200 $\mu\text{g}/\text{min}$ suggest the presence of incipient diabetic nephropathy in adult patients with type I diabetes >5 years from diagnosis.

IMPLEMENTATION AND AUDIT

Standardization of methodology and reference ranges between laboratories would be of assistance in decreasing the uncertainty associated with the various definitions of ACR and PCR.

SUGGESTIONS FOR FUTURE RESEARCH

- 1 Establish the measures of agreement between morning, random urine and 24 h urine samples.
- 2 Investigate the variability of morning, random urine and 24 h urine samples over time.
- 3 Establish the utility of gender-specific reference ranges for ACR and PCR in predicting progressive renal impairment.

REFERENCES

- Barratt A, Craig J, Cumming R. *et al. A Feasibility Study of the Early Detection and Treatment of Renal Disease by Mass Screening*. Sydney: The University of Sydney, 2003.
- Gilbert RE, Akdeniz A, Jerums G. Semi-quantitative determination of microalbuminuria by urinary dipstick. *Aust. NZ J. Med.* 1992; **22**: 334–7.
- Gerber LM, Johnston K, Alderman MH. Assessment of a new dipstick test in screening for microalbuminuria in patients with hypertension. *Am. J. Hypertens.* 1998; **11** (11 Part 1): 1321–7.
- Shaw AB, Risdon P, Lewis-Jackson JD. Protein creatinine index and Albustix in assessment of proteinuria. *BMJ. (Clin. Res. Ed)* 1983; **287**: 929–32.
- Mitchell SC, Sheldon TA, Shaw AB. Quantification of proteinuria: a re-evaluation of the protein/creatinine ratio for elderly subjects. *Age Ageing* 1993; **22**: 443–9.
- Wilson DM, Anderson RL. Protein-osmolality ratio for the quantitative assessment of proteinuria from a random urinalysis sample. *Am. J. Clin. Pathol.* 1993; **100**: 419–24.
- Ramos JG, Martins-Costa SH, Mathias MM *et al.* Urinary protein/creatinine ratio in hypertensive pregnant women. *Hypertens Pregnancy* 1999; **18**: 209–18.
- Scheid DC, McCarthy LH, Lawler FH *et al.* Screening for microalbuminuria to prevent nephropathy in patients with diabetes: systematic review of the evidence. *J. Fam Pract* 2001; **50**: 661–8.
- Verhave JC, Hillege HL, de Zeeuw D *et al.* How to measure the prevalence of microalbuminuria in relation to age and gender? *Am. J. Kidney Dis* 2002; **40**: 436–7.
- Mattix HJ, Hsu CY, Shaykevich S *et al.* Use of the albumin/creatinine ratio to detect microalbuminuria: implications of sex and race. *J. Am. Soc. Nephrol.* 2002; **13**: 1034–9.
- Jones CA, Francis ME, Eberhardt MS *et al.* Microalbuminuria in the US population. Third National Health Nutrition Examination Survey. *Am. J. Kidney Dis* 2002; **39**: 445–59.
- Garg AX, Kiberd BA, Clark WF *et al.* Albuminuria and renal insufficiency prevalence guides population screening: results from the NHANES III. *Kidney Int.* 2002; **61**: 2165–75.
- Kramer HM, Hsu CY, Curhan G. Need for sex-specific ACR. *Am. J. Kidney Dis.* 2002; **40**: 435–6.
- Ginsberg JM, Chang BS, Matarese RA *et al.* Use of single voided urine samples to estimate quantitative proteinuria. *N. Engl. J. Med.* 1983; **309**: 1543–6.
- Kristal B, Shasha SM, Labin L *et al.* Estimation of quantitative proteinuria by using the protein-creatinine ratio in random urine samples. *Am. J. Nephrol.* 1988; **8**: 198–203.
- Koopman MG, Krediet RT, Koomen GC *et al.* Circadian rhythm of proteinuria: consequences of the use of urinary protein: creatinine ratios. *Nephrol. Dial Transplant* 1989; **4**: 9–14.
- Stehouwer CD, Fischer HR, Hackeng WH *et al.* Diurnal variation in urinary protein excretion in diabetic nephropathy. *Nephrol. Dial Transplant* 1991; **6**: 238–43.
- Cundy TF, Nixon D, Berkahn L *et al.* Measuring the albumin excretion rate: agreement between methods and biological variability. *Diabet Med.* 1992; **9**: 138–43.
- Price CP, Newman DJ. *Principles and Practice of Immunoassay*, 2nd edn. London: MacMillan Reference Ltd 1997.
- Comper WD, Osicka TM, Clark M *et al.* Earlier detection of microalbuminuria in diabetic patients using a new urinary albumin assay. *Kidney Int.* 2004; **65**: 1850–5.
- Spooren PF, Lekkerkerker JF, Vermes I. Micral-Test: a qualitative dipstick test for micro-albuminuria. *Diabetes Res. Clin. Pract* 1992; **18**: 83–7.
- Marshall SM, Shearing PA, Alberti KG. Micral-test strips evaluated for screening for albuminuria. *Clin. Chem* 1992; **38**: 588–91.
- Mogensen CE, Viberti GC, Peheim E *et al.* Multicenter evaluation of the Micral-Test II test strip, an immunologic rapid test for the detection of microalbuminuria. *Diabetes Care* 1997; **20**: 1642–6.
- Parag KB, Seedat YK. The protein/creatinine index. A semiquantitative assessment of 24-hour protein excretion. *S Afr Med. J.* 1986; **69**: 42–3.
- Ralston SH, Caine N, Richards I *et al.* Screening for proteinuria in a rheumatology clinic: comparison of dipstick testing, 24 hour urine quantitative protein, and protein/creatinine ratio in random urine samples. *Ann Rheum Dis.* 1988; **47**: 759–63.
- Robert M, Sepandj F, Liston RM *et al.* Random protein-creatinine ratio for the quantitation of proteinuria in pregnancy. *Obstet Gynecol.* 1997; **90**: 893–5.
- Ruggenenti P, Gaspari F, Perna A *et al.* Cross sectional longitudinal study of spot morning urine protein: creatinine ratio, 24 hour urine protein excretion rate, glomerular filtration rate, and end stage renal failure in chronic renal disease in patients without diabetes. *BMJ.* 1998; **316**: 504–9.
- Zelmanovitz T, Gross JL, Oliveira J *et al.* Proteinuria is still useful for the screening and diagnosis of overt diabetic nephropathy. *Diabetes Care* 1998; **21**: 1076–9.
- Gatling W, Knight C, Hill RD. Screening for early diabetic nephropathy: which sample to detect microalbuminuria? *Diabet Med.* 1985; **2**: 451–5.
- Watts GF, Williams I, Morris RW *et al.* An acceptable exercise test to study microalbuminuria in type 1 diabetes. *Diabet Med.* 1989; **6**: 787–92.
- Zelmanovitz T, Gross JL, Oliveira JR *et al.* The receiver operating characteristics curve in the evaluation of a random urine specimen as a screening test for diabetic nephropathy. *Diabetes Care* 1997; **20**: 516–9.
- Ahn CW, Song YD, Kim JH *et al.* The validity of random urine specimen albumin measurement as a screening test for diabetic nephropathy. *Yonsei Med. J.* 1999; **40**: 40–5.
- Hutchison AS, O'Reilly DS, MacCuish AC. Albumin excretion rate, albumin concentration, and albumin/creatinine ratio compared for screening diabetics for slight albuminuria. *Clin. Chem* 1988; **34**: 2019–21.
- Wiegmann TB, Chonko AM, Barnard MJ *et al.* Comparison of albumin excretion rate obtained with different times of collection. *Diabetes Care* 1990; **13**: 864–71.
- Bakker AJ. Detection of microalbuminuria. Receiver operating characteristic curve analysis favors albumin-to-creatinine ratio over albumin concentration. *Diabetes Care* 1999; **22**: 307–13.
- Nathan DM, Rosenbaum C, Protasowicki VD. Single-void urine samples can be used to estimate quantitative microalbuminuria. *Diabetes Care* 1987; **10**: 414–8.
- James MA, Fotherby MD, Potter JF. Screening tests for microalbuminuria in non-diabetic elderly subjects and their relation to blood pressure. *Clin. Sci. (Lond.)* 1995; **88**: 185–90.
- Marshall SM, Alberti KG. Screening for early diabetic nephropathy. *Ann. Clin. Biochem.* 1986; **23**: 195–7.
- Eshoj O, Feldt-Rasmussen B, Larsen ML *et al.* Comparison of overnight, morning and 24-hour urine collections in the assessment of diabetic microalbuminuria. *Diabet Med.* 1987; **4**: 531–3.